

Lady Reading Hospital (LRH) Medical Teaching Institution (MTI) Peshawar

BID SOLICITATION DOCUMENTS

For

Cardiothoracic/General Surgery Disposables 2021-2022

Note: The prospective bidder is expected to examine the Bidding Documents carefully, including all Instructions, Terms & Conditions, and Specifications etc.

Failure to furnish all information required by the Bidding documents or submission of a Bid not substantially responsive to the Bidding Documents in every respect would result in the rejection of the Bid.

Contact No: 091-9211927

INTRODUCTION:

Medical Teaching Institution (MTI)/ Lady Reading Hospital Peshawar invites sealed bids from the eligible bidders for procurement of Cardiothoracic disposables for Hospital under Open Competitive Bidding under rule 6(2) (a) "Single Stage Single Envelope" bidding procedures of Khyber Pakhtunkhwa Public Procurement Regulatory Authority (KPPRA) Rules 2014.

1) INSTRUCTIONS TO BIDDERS:

- a) This Bidding procedure will be conducted in light of Khyber Pakhtunkhwa Public Procurement Regulatory Authority (KPPRA) Laws, Rules made there under, along with Standard Bidding Documents.
- b) The bids will be opened on 12 July 2021 at 10:30 AM in presence of the committee and bidders/representatives who choose to attend.
- c) Pre-bid meeting with the interested bidders will be held on 06 July 2021 at 11:00 AM in Committee room of MMD of the Institution.
- d) Any bid received after the deadline for submission of bids shall not be entertained and shall be returned unopened to the Bidder.
- e) The bid should be complete in all respect and must be signed by the bidder.
- f) All prices quoted must be in Pak Rupees (PKR) and should include all applicable taxes. If not specifically mentioned in the Quotation, it will be presumed that the prices include all the taxes. LRH Peshawar is exempted from GST.
- g) Bidders are essentially required to provide correct and latest postal/email/web addresses, phone/mobile/fax numbers for actively and timely communication.
- h) For any query, clarification regarding Services / Bid Solicitation Documents, the applicants may send a written request at least one week prior to the opening date.
- i) The Bidder may after its submission withdraw its bid prior to the expiry of the deadline prescribed for submission of bids. Withdrawn bids will be returned unopened to the Bidders.
- j) Any bid not received as per terms and conditions laid down in this document are liable to be ignored. No offer shall be considered if:
 - a. Received without earnest money;
 - b. It is received after the date and time fixed for its receipt;
 - c. The offer is ambiguous;
 - d. The offer is conditional i.e. advance payment, or currency fluctuations etc.;
 - e. The offer is from blacklisted firm in any Federal / Provincial Govt. Dept.
 - f. Only typed tender on original prescribing letter pad, sealed & signed should be submitted, the quoted Price must be preprinted and hand written quoted price will not be acceptable. The tenders must be according to hospital specification;
 - g. alternate rates (Double rates for single Items) will not be accepted.

- k) Usage of correction fluid & corrections are strictly prohibited unless duly initialed.
- 1) Any erasing / cutting etc. appearing on the offer, must be properly signed by the person signing the tender.
- m) Bids will be rejected if the Bid is in some way connected with bids submitted under names different from his own.
- n) In case of Bid Tie, the decision will be taken by making toss/draw/Recall sealed quotation in front of the bidders.
- o) Any direct or indirect effort by a bidding firm to influence this institution during the process of selection of a bidder or award of contract may besides rejection of its bid result into its disqualification from participation in the LRH Peshawar's future bids.

2) ELIGIBILITY CRITERIA:

- Bidders must give compliance to the below mentioned clauses as these are mandatory to being eligible for the bidding process. Relevant certificates must be attached.
- The bidder shall provide an undertaking that the bidder has not been declared black listed by any Governmental/ Semi-Governmental institutions.
- Bidders shall not be eligible to bid if they are under a declaration of Ineligibility for corrupt and fraudulent practices issued by any government organization in accordance with the Section 44(1) KPPRA Rules 2014.

3. GENERAL CONDITIONS: -

- 1. LRH Peshawar shall evaluate the proposal in a manner prescribed in advance, without reference to the price and reject any proposal which does not conform to the specified requirements.
- 2. Alternative bid shall not be considered and shall be rejected by the Competent Authority.
- 3. At any time prior to the deadline for submission of bids, LRH Peshawar may, for any reason, whether at its own initiative or in response to a clarification requested by a prospective Bidder, modify the bidding documents by amendment.
- 4. If a bid is not substantially responsive, it will be rejected by the Procuring Entity and may not subsequently be made responsive by the Bidder by correction of the nonconformity.
- 5. LRH MTI Peshawar may accept or reject any or all of the bids under Rule 47 of KPPRA Rules, 2014.

Lady Reading Hospital Medical Teaching institute Peshawar

Phone: 091-9211927

INVITATION FOR BIDS

Hospital Director, Medical Teaching Institute, Lady Reading Hospital Peshawar Invites sealed tenders under National Competitive Bidding for the procurement of Cardiothoracic Disposable for Hospital, under rule 6(2) (a) "single stage Single envelope procedure" of KPPRA Rules 2014, from Manufacturers / Importers / Authorized Distributors having registered Income / Sales tax, reflected on Active Taxpayer List of FBR.

The bidders are required to submit bid security **Rs-200,000**/- in the name of **Hospital Director LRH-MTI Peshawar**. Pre-bid meeting with the interested bidders will be held on 06-July 2021 at 11:00 AM at Committee room of Material Management Department of the institution.

The bids complete in all aspects must reach the undersigned by 10:00 AM on 12 July 2021 which will be opened at 10:30 AM on the same day in committee room of the Hospital in the presence of the committee and the bidders / representatives who may choose to attend.

Competent Authority reserves the right to reject any or all the bids as per provisions contained in rule 47 of KPPRA Rules 2014.

4. BID Security

Bid security Rs-200,000/- in favor of "Hospital Director LRH-MTI Hospital" The Bid Security in shape of PO (Pay Order) is not accepted.

The bid security may be forfeited:

- i) If the Bidder withdraws its bid during the period of bid validity specified by the Bidder on the Bid Form; or
- ii) In the case of a successful Bidder, if the Bidder fails to sign the contract
- iii) If a successful Bidder unable to supply the items after issuance of 3 reminders, order will be placed to the next lowest bidders at condition of risk and cost factor.

5. BID VALIDITY:

- i) The bids should be valid for 120 days.
- ii) In exceptional circumstances, LRH Hospital may solicit the Bidder's consent to an extension of the period of validity reasons shall be recorded in writing. The request and the responses there to shall be made in writing. The bid security provided shall also be suitably extended. A Bidder may refuse the request without forfeiting its bid security. A Bidder granting the request will not be required nor permitted to modify its bid, except as provided in the bidding document.

CARDIOTHORACIC DISPOSABLES LIST FOR THE YEAR 2021-22.

S.NO	Product
1	New Titanium Linear Cutter 55mm
2	New Titanium Linear Cutter 75mm
3	New Titanium Linear Cutter Reload 55mm
4	New Titanium Linear Cutter Reload 75mm
5	Enseal G2 Curved Articulating Tissue Sealer 5mm D,35cm
6	GEN11.Harmonic Connector
7	Echelon Flex 45mm Standard Articulating Stapler
8	Echelon Flex 60mm Standard Articulating Stapler
9	Echelon Endopath Standard Reload White 45mm
10	Echelon Endopath Standard Reload Green 45mm
11	Echelon Endopath Standard Reload Blue 45mm
12	Echelon Endopath Standard Reload White 60mm
13	Echelon Endopath Standard Reload Blue 60mm
14	Echelon Endopath Standard Reload Green 60mm
15	Harmonic Blue Handpiece
16	Handpiece Harmonic Scalpel
17	Harmonic Focus Shear Curved 9cm
18	Enseal G2 Curved Tissue Sealer 5mm D,35cm Shaft Length
19	Harmonic ACE+ Shear 5mm D,23CM Length
20	Harmonic ACE+ Shear 5mm D,36CM Length

S	.No	Nomenclature	
1		Thoracoports with sized 5.5mm, 10.5mm, 11.5mm, and 15mm	
2		Laparoscopic/VATS surgery Liner cutter with facility to reload 30mm, 45mm, 60mm Straight /Roticulator Cartridges in one gun without integrated knife with 360* rotation & 0-45* articulation in both directions	

The EEA'm circular stapler with DST Series' technology embodies our commitment to innovation and collaboration with the surgical community. Our circular staplers incorporate more than three decades of experience and input from leading bariatric, colorectal, and general surgeons throughout the world. DST Series' EEA's staplers incorporate proprietary directional stapling technology (DST), 3.5 mm staples, and a host of other unique features. The DST Series' EEA' OrVil' 21 mm device is a DST Series' EEA' 21 mm anvil assembly, and the DST Series' EEA' OrVil' 25 mm device is a DST Series' EEA' 25 mm anvil assembly, both with the anvil head secured in the tilted position. The anvil assembly is mounted on a 90 cm long PVC delivery tube and is secured to the tube with a suture. Only for use with the DST Series' EEA' XL 21 or 25 mm stapler. Tri-Staple' technology reloads (Tan, Purple, Black) are designed to work in harmony with the natural properties of tissue to optimize performance before, during and after stapling. With its stepped cartridge face, Tri-Staple' technology delivers graduated compression, optimizing interaction between the tissue and stapler because it Generates less stress on tissue during compression and clamping1,2,1,‡ Allows greater perfusion into the staple line3,1,‡ Provides superior performance in variable thicknesses2,4,5,6,7,† Tri-Staple' technology improves staple line strength, leak resistance and hemostasis when compared to the universal reloads. Disposable 5mm Endoscopic grasper with Ratcheted Mechanism Disposable 5mm Endoscopic daryland Dissector Disposable 5mm Endoscopic daryland Dissector Roticulator Endoscopic Dissector Sim with On-00ff ratchet switch and 360* rotation Roticulator Endoscopic Dissector Sim with On-00ff ratchet switch and 360* rotation Roticulator Endoscopic Dissector Sim with On-00ff ratchet switch and 360* rotation Roticulator Endoscopic Dissector Sim with On-00ff ratchet switch and 360* rotation Roticulator Endoscopic Grasper Sim with On-00ff ratchet switch a		
The DST Series™ EEA™ OrVil™ 21 mm device is a DST Series™ EEA™ 25 mm anvil assembly, and the DST Series™ EEA™ OrVil™ 25 mm device is a DST Series™ EEA™ 25 mm anvil assembly, both with the anvil head secured in the tilted position. The anvil assembly is mounted on a 90 cm long PVC delivery tube and is secured to the tube with a suture. Only for use with the DST Series™ EEA™ XL 21 or 25 mm stapler. Tri-Staple™ technology reloads (Tan,Purple,Black) are designed to work in harmony with the natural properties of tissue to optimize performance before, during and after stapling. With its stepped cartridge face, Tri-Staple™ technology delivers graduated compression, optimizing interaction between the tissue and stapler because it. Generates less stress on tissue during compression and clamping1,2,1,‡ Allows greater perfusion into the staple line3,1,‡ Provides superior performance in variable thicknesses2,4,5,6,7,† Tri-Staple™ technology improves staple line strength, leak resistance and hemostasis when compared to the universal reloads. Disposable 5mm Endoscopic Maryland Dissector Disposable 5mm Endoscopic Maryland Dissector Disposable 5mm Endoscopic durved scissor like a Metzenbaum with cautry Endo lung 12mm Tringale shape designed for grasping lung and other tissues Defined Clinch 5mm atraumatic grasper with On/Off ratchet switch and 360* rotation Roticulator Endoscopic Dissector 5mm with 0-80 Degree articulating tip and 360* rotation Roticulator Endoscopic Grasper 5mm with 0-80 Degree articulating tip and 360* rotation Roticulator Endoscopic Grasper 5mm with 0-80 Degree articulating tip and 360* rotation Disposable suture loading units for Endostitch instrument: polyester, silk, nylon, absorbable suture Endo Catch Specimen Removal Kit: 10mm with leak proof and impervious material to cancer cells of 0.5 microns/ pretied purse string on pouch. Endo Raga specimen retrieval system sizes:3x6 & 5x8 Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin plac	4	innovation and collaboration with the surgical community. Our circular staplers incorporate more than three decades of experience and input from leading bariatric, colorectal, and general surgeons throughout the world. DST Series™ EEA™ staplers incorporate proprietary directional stapling technology (DST),
natural properties of tissue to optimize performance before, during and after stapling. With its stepped cartridge face, Tri-Staple™ technology delivers graduated compression, optimizing interaction between the tissue and stapler because it Generates less stress on tissue during compression and clamping1,2,†,‡ Allows greater perfusion into the staple line3,1,‡ Provides superior performance in variable thicknesses2,4,5,6,7,† Tri-Staple™ technology improves staple line strength, leak resistance and hemostasis when compared to the universal reloads. 6 Disposable 5mm Endoscopic grasper with Ratcheted Mechanism 7 Disposable 5mm Endoscopic Maryland Dissector 8 Disposable 5mm Endoscopic Maryland Dissector 9 Endo lung 12mm Tringale shape designed for grasping lung and other tissues 10 Endo Clinch 5mm atraumatic grasper with On/Off ratchet switch and 360* rotation 11 Roticulator Endoscopic Dissector 5mm with 0-80 Degree articulating tip and 360* rotation 12 Roticulator Endoscopic Grasper 5mm with 0-80 Degree articulating tip and 360* rotation 13 Endostitch Instrument 10 mm for use with double point needle suture combination 14 Disposable suture loading units for Endostitch instrument: polyester, silk, nylon, absorbable suture 15 Endo Catch Specimen Removal Kit: 10mm with leak proof and impervious material to cancer cells of 0.5 microns/ pretied purse string on pouch. 16 Endo Retract 10mm three finger Retracter 17 Endo Bag specimen retrieval system sizes:3x6 & 5x8 Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 30mm. Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm Liner Stapler with single lever operation & thumb push button for op	4	The DST Series™ EEA™ OrVil™ 21 mm device is a DST Series™ EEA™ 21 mm anvil assembly, and the DST Series™ EEA™ OrVil™ 25 mm device is a DST Series™ EEA™ 25 mm anvil assembly, both with the anvil head secured in the tilted position. The anvil assembly is mounted on a 90 cm long PVC delivery tube and is secured to the tube with a suture.
Disposable 5mm Endoscopic Maryland Dissector Disposable 5mm Endoscopic curved scissor like a Metzenbaum with cautry Endo lung 12mm Tringale shape designed for grasping lung and other tissues Endo Clinch 5mm atraumatic grasper with On/Off ratchet switch and 360* rotation Roticulator Endoscopic Dissector 5mm with 0-80 Degree articulating tip and 360* rotation Roticulator Endoscopic Grasper 5mm with 0-80 Degree articulating tip and 360* rotation Endostitch Instrument 10 mm for use with double point needle suture combination Disposable suture loading units for Endostitch instrument: polyester, silk, nylon, absorbable suture Endo Catch Specimen Removal Kit: 10mm with leak proof and impervious material to cancer cells of 0.5 microns/ pretied purse string on pouch. Endo Retract 10mm three finger Retracter Endo Bag specimen retrieval system sizes:3x6 & 5x8 Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 30mm. Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm	5	natural properties of tissue to optimize performance before, during and after stapling. With its stepped cartridge face, Tri-Staple™ technology delivers graduated compression, optimizing interaction between the tissue and stapler because it: Generates less stress on tissue during compression and clamping1,2,†,‡ Allows greater perfusion into the staple line3,†,‡ Provides superior performance in variable thicknesses2,4,5,6,7,† Tri-Staple™ technology improves staple line strength, leak resistance and hemostasis when
Disposable 5mm Endoscopic curved scissor like a Metzenbaum with cautry Endo lung 12mm Tringale shape designed for grasping lung and other tissues Endo Clinch 5mm atraumatic grasper with On/Off ratchet switch and 360* rotation Roticulator Endoscopic Dissector 5mm with 0-80 Degree articulating tip and 360* rotation Roticulator Endoscopic Grasper 5mm with 0-80 Degree articulating tip and 360* rotation Endostitch Instrument 10 mm for use with double point needle suture combination Disposable suture loading units for Endostitch instrument: polyester, silk, nylon, absorbable suture Endo Catch Specimen Removal Kit: 10mm with leak proof and impervious material to cancer cells of 0.5 microns/ pretied purse string on pouch. Endo Retract 10mm three finger Retracter Endo Bag specimen retrieval system sizes:3x6 & 5x8 Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 30mm. Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm	6	
9 Endo lung 12mm Tringale shape designed for grasping lung and other tissues 10 Endo Clinch 5mm atraumatic grasper with On/Off ratchet switch and 360* rotation 11 Roticulator Endoscopic Dissector 5mm with 0-80 Degree articulating tip and 360* rotation 12 Roticulator Endoscopic Grasper 5mm with 0-80 Degree articulating tip and 360* rotation 13 Endostitch Instrument 10 mm for use with double point needle suture combination 14 Disposable suture loading units for Endostitch instrument: polyester, silk, nylon, absorbable suture 15 Endo Catch Specimen Removal Kit: 10mm with leak proof and impervious material to cancer cells of 0.5 microns/ pretied purse string on pouch. 16 Endo Retract 10mm three finger Retracter 17 Endo Bag specimen retrieval system sizes:3x6 & 5x8 Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 30mm. Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm	7	Disposable 5mm Endoscopic Maryland Dissector
Endo Clinch 5mm atraumatic grasper with On/Off ratchet switch and 360* rotation Roticulator Endoscopic Dissector 5mm with 0-80 Degree articulating tip and 360* rotation Roticulator Endoscopic Grasper 5mm with 0-80 Degree articulating tip and 360* rotation Endostitch Instrument 10 mm for use with double point needle suture combination Disposable suture loading units for Endostitch instrument: polyester, silk, nylon, absorbable suture Endo Catch Specimen Removal Kit: 10mm with leak proof and impervious material to cancer cells of 0.5 microns/ pretied purse string on pouch. Endo Retract 10mm three finger Retracter Endo Bag specimen retrieval system sizes:3x6 & 5x8 Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 30mm. Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology	8	Disposable 5mm Endoscopic curved scissor like a Metzenbaum with cautry
11 Roticulator Endoscopic Dissector 5mm with 0-80 Degree articulating tip and 360* rotation 12 Roticulator Endoscopic Grasper 5mm with 0-80 Degree articulating tip and 360* rotation 13 Endostitch Instrument 10 mm for use with double point needle suture combination 14 Disposable suture loading units for Endostitch instrument: polyester, silk, nylon, absorbable suture 15 Endo Catch Specimen Removal Kit: 10mm with leak proof and impervious material to cancer cells of 0.5 microns/ pretied purse string on pouch. 16 Endo Retract 10mm three finger Retracter 17 Endo Bag specimen retrieval system sizes:3x6 & 5x8 Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 30mm. Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology	9	Endo lung 12mm Tringale shape designed for grasping lung and other tissues
12 Roticulator Endoscopic Grasper 5mm with 0-80 Degree articulating tip and 360* rotation 13 Endostitch Instrument 10 mm for use with double point needle suture combination 14 Disposable suture loading units for Endostitch instrument: polyester, silk, nylon, absorbable suture 15 Endo Catch Specimen Removal Kit: 10mm with leak proof and impervious material to cancer cells of 0.5 microns/ pretied purse string on pouch. 16 Endo Retract 10mm three finger Retracter 17 Endo Bag specimen retrieval system sizes:3x6 & 5x8 Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 30mm. Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology	10	Endo Clinch 5mm atraumatic grasper with On/Off ratchet switch and 360* rotation
13 Endostitch Instrument 10 mm for use with double point needle suture combination 14 Disposable suture loading units for Endostitch instrument: polyester, silk, nylon, absorbable suture 15 Endo Catch Specimen Removal Kit: 10mm with leak proof and impervious material to cancer cells of 0.5 microns/ pretied purse string on pouch. 16 Endo Retract 10mm three finger Retracter 17 Endo Bag specimen retrieval system sizes:3x6 & 5x8 Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 30mm. Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology	11	Roticulator Endoscopic Dissector 5mm with 0-80 Degree articulating tip and 360* rotation
Disposable suture loading units for Endostitch instrument: polyester, silk, nylon, absorbable suture Endo Catch Specimen Removal Kit: 10mm with leak proof and impervious material to cancer cells of 0.5 microns/ pretied purse string on pouch. Endo Retract 10mm three finger Retracter Endo Bag specimen retrieval system sizes:3x6 & 5x8 Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 30mm. Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology Isomorphic for pening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology	12	Roticulator Endoscopic Grasper 5mm with 0-80 Degree articulating tip and 360* rotation
Endo Catch Specimen Removal Kit: 10mm with leak proof and impervious material to cancer cells of 0.5 microns/ pretied purse string on pouch. Endo Retract 10mm three finger Retracter Endo Bag specimen retrieval system sizes:3x6 & 5x8 Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 30mm. Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology automatic / manual pin placement slim anvil, longer shaft, Directional staple technology	13	Endostitch Instrument 10 mm for use with double point needle suture combination
cancer cells of 0.5 microns/ pretied purse string on pouch. Endo Retract 10mm three finger Retracter Endo Bag specimen retrieval system sizes:3x6 & 5x8 Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 30mm. Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology	14	
Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 30mm. Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology	15	
Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 30mm. Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology	16	Endo Retract 10mm three finger Retracter
automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 30mm. Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology	17	Endo Bag specimen retrieval system sizes:3x6 & 5x8
automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 45mm Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology	18a	automatic / manual pin placement slim anvil, longer shaft, Directional staple technology
18c automatic / manual pin placement slim anvil, longer shaft, Directional staple technology	18b	automatic / manual pin placement slim anvil, longer shaft, Directional staple technology
	18c	Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology

18d	Liner Stapler with single lever operation & thumb push button for opening the jaws, automatic / manual pin placement slim anvil, longer shaft, Directional staple technology incorporating flat wire staples for better control in sizes of 90mm
19a	Cartridges for liner Stapler with directional staple technology incorporating flat wire staples in sizes of 30mm, with white 2.5/ blue 3.5& green 4.8 leg lengths
19b	Cartridges for liner Stapler with directional staple technology incorporating flat wire staples in sizes of 45mm, with white 2.5/ blue 3.5& green 4.8 leg lengths
19c	Cartridges for liner Stapler with directional staple technology incorporating flat wire staples in sizes of 60mm with white 2.5/ blue 3.5& green 4.8 leg lengths
19d	Cartridges for liner Stapler with directional staple technology incorporating flat wire staples in sizes of 90mm with white 2.5/ blue 3.5& green 4.8 leg lengths
20a	Liner Cutter with Adjustable lever operation from both hand side, with rectangular wire staples/flat wire staples, gun without integrated knife, in sizes of 60mm.
20b	Liner Cutter with Adjustable lever operation from both hand side, with rectangular wire staples/flat wire staples, gun without integrated knife, in sizes of 80mm.
20c	Liner Cutter with Adjustable lever operation from both hand side, with rectangular wire staples/flat wire staples, gun without integrated knife, in sizes of 100mm
21a	Cartridges for liner cutter, incorporating flat wire staples in sizes of 60mm with white 2.5mm, blue 3.8mm& green 4.8mm leg lengths, with integrated knife in each cartridge/reload
21b	Cartridges for liner cutter, incorporating flat wire staples in sizes of 80mm with white 2.5mm, blue 3.8mm& green 4.8mm leg lengths, with integrated knife in each cartridge/reload
21c	Cartridges for liner cutter, incorporating flat wire staples in sizes of 100mm with white 2.5mm, blue 3.8mm& green 4.8mm leg lengths, with integrated knife in each cartridge/reload
22	LigaSure Maryland jaw 37 cm Sealer/Divider Laparoscopic
23	LigaSure Maryland jaw 23 cm Sealer/Divider open
24	LigaSure Blunt Tip 5mm Laparoscopic
25	LigaSure Atlas 10mm 20cm
26	LigaSure Atlas 10mm 37 cm
27	LigaSure small Jaw for Thyroidectomy.
28	Liga sure Impact Open Instrument
29	Chest Bottle
30	Bronchile Tube Right & Left all sizes

7. Payment:

- a. No advance payment will be permissible.
- b. The payment will be made after successful supply and inspection.